

LOC255167-hypothetical LOC255167 Gene

Reference: PTXBC106918

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC255167
GENE ID	255167
ORIGIN SPECIES	Human
ACCESSIONS	BC106918
VECTOR	pDONR223
NCBI SYMBOL	LOC255167
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggacctcgacaacagcgcttgcggtggcccttctgcagctctccacccgcgcgtgctcgggacactgggcttgaaggcgccagttccagccgcacgtgcagcctaggaggcgccagggcgctgtccctgtggagggaagggccggcaacc
tgcaaaagacgtcacagaaacttgcggcttgcggccgatgcgcgtcaggtgctccctccaggtggcagaggagggaggggtgcgcctggtggggcctgagagccaccgggctgcgtgcgagtggtggacaccggccgcttctgctctacagcct
gggtctgggaccttggaaaccgcttggacgaagctga