



Reference: PTXBC106936

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CNCA2; CNG-4; CNG4; CNG5; CNGB2; OCNC2; OCNCBETA; OCNCb; cyclic nucleotide-gated cation channel alpha-4; CNG channel alpha-4; cyclic nucleotide gated channel beta 2; cyclic nucleotide gated channel alpha 4		

For research use only.

GENE DESCRIPTION	cyclic nucleotide gated channel alpha 4
GENE ID	1262
ORIGIN SPECIES	Human
ACCESSIONS	BC106936
VECTOR	pDONR223
NCBI SYMBOL	CNGA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr