

BEX5-brain expressed, X-linked 5 Gene

Reference: PTXBC106955

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein BEX5; NGFRAP1L1; BEX family member 5; NGFRAP1-like 1; NGFRAP1-like protein 1; brain-expressed X-linked protein 5; nerve growth factor receptor-associated protein 2; brain expressed X-linked 5		

For research use only.

GENE DESCRIPTION	brain expressed, X-linked 5
GENE ID	340542
ORIGIN SPECIES	Human
ACCESSIONS	BC106955
VECTOR	pDONR223
NCBI SYMBOL	BEX5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaatgtcccaaggaaaacaaagtgtggagaaggcccgagtcagaaatgaagcccccgttaggaggtggtaataccaggagcctggaggaaatgttaaaggggtttgggtccacgtgccccgggtttggagaggatgtgcccaataggcttgataacattgatatg
atagatggagatggagatgatatggaacgggttcagaggagatgagagagctaaaggagaaatagggaacttcagttgaggtacagtcgtcgatcttataggggaccctcctcaccatgatcatcatgatgagtttgccttatgccttga