

PRDM7-PR domain containing 7 Gene

Reference: PTXBC107033

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PFM4; ZNF910; PR domain 7; PR domain containing 7; PR domain zinc finger protein 7; PR-domain family protein 4; PR/SET domain 7		

For research use only.

GENE DESCRIPTION	PR domain containing 7
GENE ID	11105
ORIGIN SPECIES	Human
ACCESSIONS	BC107033
VECTOR	pDONR223
NCBI SYMBOL	PRDM7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtcagaactttctcattgacagctgtgctcatgggccccctacattgttaaaggacagtgcagtggaacaaggcgatccaacgttcagccctcagctgtccccggggctgagaattgggcatcaggcatccctcaggctgggcttgagatggaacgagggcatctgatctgcccactgggtctgcactttggccctatgagggccgaattacagaagacgaaggcagcccaacagtggatattcctggctaatcaccaaggggagaaactgctatgagatgtggatggaaaagataaatcctcgcccaactggatgagaaccaagccagagatccatccatgtccctcatgctgtctggccttttcaagtaaaaaatttctcagtcacatgtggaacgaatcactcctctcagaacttccaggaccatctgcaagaaaacttctccaaccagagaatccctgcccagggatcagaatcaggagcggcaatattctga