

BAGE-B melanoma antigen Gene

Reference: PTXBC107038

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BAGE1; CT2.1; B melanoma antigen 1; antigen MZ2-BA; cancer/testis antigen 2.1; cancer/testis antigen family 2, member 1; B melanoma antigen		

For research use only.

GENE DESCRIPTION	B melanoma antigen
GENE ID	574
ORIGIN SPECIES	Human
ACCESSIONS	BC107038
VECTOR	pDONR223
NCBI SYMBOL	BAGE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccagagcgggttttctggcattgtctgccagctgtccaagccaggctgatgaaggaggagtcctctgtgtgagctggagggtggagcctgaagacggcacagctctgtgttcattcttga