

## IL1F6-interleukin 1 family, member 6 (epsilon) Gene

Reference: PTXBC107043

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL1F6; FIL1; FIL1(EPSILON); FIL1E; IL-1F6; IL1(EPSILON); interleukin-36 alpha; FIL1 epsilon; IL-1 epsilon; IL-1F6 (FIL-1-epsilon); interleukin 1 family, member 6 (epsilon); interleukin-1 epsilon; interleukin-1 family member 6; interleukin 36, alpha		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	interleukin 1 family, member 6 (epsilon)
<b>GENE ID</b>	27179
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC107043
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	IL1F6
<b>LENGTH</b>	Proteogenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atggaaaagcattgaaaatgacacacctcagcaggggagcattcaggatcaatcatcgggtgtggttctcaggaccagcgtcatagcagtcgccgaggaagaccgtatgtctccagtcactattgccttaatctcatgccgacatgtggagacccttgagaaaagacagaggggaacccatctacctggcctgaatggactcaatctctgcctgatgtgtcctaaagtctgggacagccacactgcagctgaaggaaggaatataatggattgtacaaccaaccgagcctgtgaagtccttctctctaccacagccagagtgccaggaactccacctcgagctctggtcttccctggctgctcagctctgaaggagcctgctctcctcctaccgaactggggaagccaactactgacttgggttaactatgctgttttaa