

PATE-expressed in prostate and testis Gene

Reference: PTXBC107044

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PATE; prostate and testis expressed protein 1; expressed in prostate and testis; prostate and testis expressed 1		

For research use only.

GENE DESCRIPTION	expressed in prostate and testis
GENE ID	160065
ORIGIN SPECIES	Human
ACCESSIONS	BC107044
VECTOR	pDONR223
NCBI SYMBOL	PATE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacaagtcctcttctggaactccccatcctgctctgctctttagggcattatctggatcactttcaatgagaaatgatgcagtgatagaaattgtcgggtgtaggatgtgccacctccagttcccaggagaaaagtgtccagaggagagggaatgcacagcaacaacagaagag
gcctgcagtggttgaaggatgttcaaaaggatgtaaccttcatgggctgcctaaagaactgtgctgatgtgaaaggcataagggtggagtgctatttggtaacttcagggtgctgcaggagccatgacctgtgcaatgaagaccttttag