

DKK4-dickkopf homolog 4 (*Xenopus laevis*) Gene

Reference: PTXBC107047

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DKK-4; dickkopf-related protein 4; dickkopf-4; hDkk-4; dickkopf WNT signaling pathway inhibitor 4		

For research use only.

GENE DESCRIPTION	dickkopf homolog 4 (<i>Xenopus laevis</i>)
GENE ID	27121
ORIGIN SPECIES	Human
ACCESSIONS	BC107047
VECTOR	pDONR223
NCBI SYMBOL	DKK4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggcgccgctgctgctgggctgagctgctgctccccctgggagctgctgctggaactcaacaacatcaggagctgctgacctgcatggggcccggaagggtcacagtgcctgtctgacacggactgcaataccagaaagtctgctccagccccgcgatgagaagccgttctgtgctacatgctgctgggtgctggaggaggtgccagcgagatgccatgctgctgctgggacactctgtgtgaacgatgtttgtactacgatggaagatgcaaccccaatattagaaggcagcttgatgagcaagatggcacacatgcagaaggacaactgggcacccagttccaggaaaaccaacccaaaaggagccaagtattaaagaatcacaaggcagggaagggaagaagggtgtctgagaacttttactgtggccctggacttttctgtgctgctgctcattttggacgaaaatttgaagccagtccttttggaggacaggtctgctccagaagagggcataaagacactgctcaagctcagaaatcttcagcgttgcgactgtggccctggactactgtgtcgaagccaattgaccagcaatcggcagcatgctcgattaagagtagtgcacaaaatagaaaagctataa