

R3HDML-R3H domain containing-like Gene

Reference: PTXBC107048

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	cysteine-rich secretory protein R3HDML; peptidase inhibitor R3HDML; dJ881L22.3; R3H domain containing like		

For research use only.

GENE DESCRIPTION	R3H domain containing-like
GENE ID	140902
ORIGIN SPECIES	Human
ACCESSIONS	BC107048
VECTOR	pDONR223
NCBI SYMBOL	R3HDML
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctctgccagcaccgtggcctggcagcctgtcttctggcctggcagcagtgaaacgcttgataatgctaaccagccccggccagccccgagagcacgctatgcggctcctgagtgccctggaggtgccaggtaccgccggaagcggccatctctg
tgagagacatgaatgccttactggattatcacaccacatccggccaggtgtgtaccacctgccgccaacatggaatacatggtctgggacaagcgtgcccagggctgccgaagcctggccaccagtgcatctgggcacatgggcttcacagctgatgagatactggccag
aacctctccatccattctggccaggtaccggtccgtagtggtatctcatgaagtcctggtctgaggagaagtggcattactgtttccggcccaagggaactgtaaccacactgccctggcgtgcgatggccacactgtccattataccagatggtgtgggcatctccaatcggtg
ggctgtgccatccacactgtatgacatcagtgctggggcaacacctggcatcggcgccatacctggtctgcaactatgccattaagggaactggattggcgagtcctccgtacaagatgggaaagccgtgctcctctgtccccccagttatcaaggcagctgcaatagcaacatgt
gcttcaaggggctgaatccaacaagttcacgtggttctga