

FKSG83-FKSG83 Gene

Reference: PTXBC107050

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	FKSG83
GENE ID	83954
ORIGIN SPECIES	Human
ACCESSIONS	BC107050
VECTOR	pDONR223
NCBI SYMBOL	FKSG83
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgggtcattataagcaccgtagccaccagcaaatccactaatgttgacttgatattctctgatatactgtaaaacaagactgtttatatactaccatgcaatgataattcaaaagtccatgtttactcaaggtttctctctgttctctcatgatttgactagtttctacatgggtgactattttcttaa
ggcatcaaaagacggctcttcatattcacagatcagctgtctccacaatctttctgaaacagaagccagtcagttttctctctgtgtgagccgctgtttttctttatgggactaacgtgtgtcttccacttatattagctccatgatgagaacagccctaatcctggagaacattactagtttg
tttcatcttctaccacaattatctgtactttgatattgatttaattatgataataaagtgtccagattaacctgtgccattttaaatatgattatttccctgcttaccacacagcttctattgcatcaaaagaatcctgataaaaattatfttaaccacaagagaatgtgttttaa