

VN1R3-vomer nasal 1 receptor 3 Gene

Reference: PTXBC107074

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	vomer nasal 1 receptor 3
GENE ID	317702
ORIGIN SPECIES	Human
ACCESSIONS	BC107074
VECTOR	pDONR223
NCBI SYMBOL	VN1R3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctccaaggatttgcgaataggaatgatcttattatctcagatcatggcggatctctgggaattcttctctaccactatagttctctgttccacagaggtatgttacagtcacagatctgattctcaagcactgaccatagccaactccttggttatactctaaaggaatccac
aaacaatggctgcttttgggtgaaagattccctcagtgataatggatgcaaatgtgtttatgttcacagagtgaggcaggctgtgtgcaactgagtgcttccagggtcatcaccatcagccccagtgga