

NTF3-neurotrophin 3 Gene

Reference: PTXBC107075

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HDNF; NGF-2; NGF2; NT-3; NT3; neurotrophin-3; nerve growth factor 2; neurotrophic factor		

For research use only.

GENE DESCRIPTION	neurotrophin 3
GENE ID	4908
ORIGIN SPECIES	Human
ACCESSIONS	BC107075
VECTOR	pDONR223
NCBI SYMBOL	NTF3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccatctgttttatgtgatatttctcgttctcgttgcacccaaggttaacaacatggatcaaaaggagtttgcagaagactcgtcaattccctcattatgaagctgatccaggcagatatttgaataaagctctccaagcagatggtggacgttaaggaaattaccagagcacct
gccccaaagctgaggtcctcccgagagccggagcggggaggcccgccaagtcagcattccagccagtgaattgcaatggacaccgaactgctgcgacaacagagacgtacaactcaccgctgctgagcgacagcacccttggagccccgccttgtatctcatggagg
attacgtgggcagccccgtggtggcgaacagaacatcacggcgaaacgggtacgcgagcataagagtcaccgagggaggtactcgtatgtgacagtgagagctctgtgggtgaccgacaagtcacggccatcgacattcggggacaccaggtcacggtgctggggagatcaaa
acgggcaactctcctgtcaacaataattttatgaaacgcgatgtaaggaaagccaggccggtcaaaaacggttcgaggggtattgatataaacactggaactctcagtcgaaaacatcccaaacctacgtccgagcactgacttcagagaacaataaactgtgggtgctgcgggtgatac
ggatagacacgtcctgtgtgtgtccttgcgagaaaaatcggaagaacatga