

APOOL-apolipoprotein O-like Gene

Reference: PTXBC107096

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CXorf33; FAM121A; Mic27; UNQ8193; MICOS complex subunit MIC27; AAIR8193; family with sequence similarity 121A; apolipoprotein O like		

For research use only.

GENE DESCRIPTION	apolipoprotein O-like
GENE ID	139322
ORIGIN SPECIES	Human
ACCESSIONS	BC107096
VECTOR	pDONR223
NCBI SYMBOL	APOOL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccatcaggatgggaaactgacaaccatgcctgcaggtctgatatgcatctgtaagtgtacatgcagccaaacaggaatccaaaagcagctagtgaaccagagcagctcccatatatactgcaccaccgctccagctctaaatattgtaagagcagcctggctatt
acaaatgggctttgctccatccgactgcaactggtgttacattggctggtgcaagggtgtttatgtctttgtaaaaaatgggataatggatacagtaacaattggaaaagatgcttatgtctatctgaagaatctctcgcagattttctccgaaaatgggagttattacagtttcaggattggcg
ggcttggttcagcgagaaaagggtccaagttaagaaaattactatctctgggactggccactttaggagcaactgtttgctacccagttcagtcctgaataattgctaaggtaacagcaaaaaagggtatgctacaagccagcaaatgttgagcagttaaatcattgtggacaaaaagca
gcaaaagagagtcactccctaacctaagaaaaaactaagctaggtatctctccgaaatagaagtacgtgcaaaaacaactcacgtcttgaacactcagtgcccttgccaacagaaactcagctctgaagcaaaagaccaaatcagaatccacttcaggtgcaaccagtttatgcctgac
cccaagctcatggtacacggcagtcacccccaagaagatagatatgtacagcactagaagctga