

UPK1A-uroplakin 1A Gene

Reference: PTXBC107098

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TSPAN21; UP1A; UPIA; UPKA; uroplakin-1a; tetraspanin-21; tspan-21; uroplakin Ia; uroplakin 1A		

For research use only.

GENE DESCRIPTION	uroplakin 1A
GENE ID	11045
ORIGIN SPECIES	Human
ACCESSIONS	BC107098
VECTOR	pDONR223
NCBI SYMBOL	UPK1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtctgcggcagcagcggaggccgagaaggatctccagttgtgtggcctgctagttgtgggcaatattctgtctgagccctgtccctgtttgtctgagaccatagggtgacagccgaccagtaccgtgtataccactgatggagctcaggaaggatgacgtcttc
gctgtgctggtgattgccatctctgcggcttctctcttctcatggtagccagttttgtgtgtgggtgcccgaactctgcccggccggctccatggtctcactgacctgtgtctcatgctcactcttcgagtgccgctcctgcatcagctctacacccacgtgactacatgggttccaa
cccatccctgataccaagcagatgctgacctctacagcgggacaccgaccaggccaggagctgacccgctctgggaccgctcatgattgagcaagaatgctgtggcacatctgtcccatggactgggtgaacttcacgtcagccttcggggccgaccctccggaggtgtgt
tccccggccccactgtgctgtcggcgacgggaaacttcacccctcaacgaggagggtgcccgtgggcacatggactactgttccaagggtgcttcgaacacatggccacgccatcgacagctacagtgagggtatctctgtgtttgggtttgccatcctgatgtggac
gtccccgtcatgctgatagccatgtatttctacaccagctctga