

C21orf84-chromosome 21 open reading frame 84 Gene

Reference: PTXBC107109

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf84; NCRNA00313; long intergenic non-protein coding RNA 313		

For research use only.

GENE DESCRIPTION	chromosome 21 open reading frame 84
GENE ID	114038
ORIGIN SPECIES	Human
ACCESSIONS	BC107109
VECTOR	pDONR223
NCBI SYMBOL	C21orf84
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacaactctgagcaagtcaaaacaacctcagcagcgggaaacctcgatgaacaacaccggagaggtttctcgatcctctttgtgacttcgagtgaagtgtggaagcctggaagcacttagaccctgcctgtgccgcagaccctgggagcaccacacccccgatcga
gaagcccgatgggaaagaatgtctgggaaggacgtcatgcccgaggcttgcgtga