

NSUN5B-NOL1/NOP2/Sun domain family, member 5B Gene

Reference: PTXBC107112

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	NOL1/NOP2/Sun domain family, member 5B
GENE ID	155400
ORIGIN SPECIES	Human
ACCESSIONS	BC107112
VECTOR	pDONR223
NCBI SYMBOL	NSUN5B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccacgctgctggcctgggtggcgtctcctgctgagctgggtgaggaggacttctggcgtctcccccttagatccgcgtatcgtgaggtccactatgtcctgctggatccttctgcagtggctcgggtatgccgagcagacagctggaggatccggggcaggacaccta
gcccgggtgctgctgcatgccctggcagggttcagcagcagccctgtgccacgcgtcacttccctccctgcagcggtcgtctactccatgtgctcccttgcaggaggagaatgaagacatggtaccagatgcgtgcagcagaacccgggccttcaggtactccgccc
tgctgccccggccccaccgaggcctgagcacgttccccgggtgccgagcactgctccgggttcccccaagaccacgcttagcgggtggtcttctgtgctgtaattgaacgggtcagatgccgacgtga