

WDR86-WD repeat domain 86 Gene

Reference: PTXBC107125

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WD repeat-containing protein 86; WD repeat domain 86		
For research use only.			

GENE DESCRIPTION	WD repeat domain 86
GENE ID	349136
ORIGIN SPECIES	Human
ACCESSIONS	BC107125
VECTOR	pDONR223
NCBI SYMBOL	WDR86
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgggagttccggggccaccgaactgcgtgctgaccctagcctactgtgcccgtgggacgtcccccagcactccctgcgcggaggaggccggcgggggggcttctgtgacgggcagcacagatggcacagccaagggtggcagggtggccagcggtgctgccacc
agacgtgcggggccacacgggtgcagtgctgtgcttagtgctagacacggccgacacggccttcacaggcagcaccgacgccaccatccgtgcttgggacatcctgagtgaggagcagctgcgggtgttccgggagcaccggggctccgtcatctgtctggagctggtgaac
cgactcgtgactctggcagcgggacaggaccgtcaagtgtggtgtggcagacacaggggagtggtgtgcgcacgttcacggccacagacgaacgtgagcggccctcaagtaccacggggacacctgttcacgggcagcggggacgcttgcggccgggcttcgacgcgcagt
ctggagagctgcggagggtgttccggggccacattcatcatcaactgcatccaggtgcacggccaggtgctctacacggcctgcacgacggcgccctgcgcctctgggacgtgcggggctccgaggtgccccgggccccctcgcccatgagcagcctctcgcgctcttca
gcaacaagggtgggtgcggccggccctgcagccggcctga