



Reference: PTXBC107127

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SAS1B; astacin-like metalloendopeptidase; astacin-like metallo-endopeptidase (M12 family); astacin-like metalloendopeptidase (M12 family); oocyte astacin; sperm acrosomal SLLP1 binding; astacin like metalloendopeptidase		

For research use only.

GENE DESCRIPTION	astacin-like metallo-endopeptidase (M12 family)
GENE ID	431705
ORIGIN SPECIES	Human
ACCESSIONS	BC107127
VECTOR	pDONR223
NCBI SYMBOL	ASTL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]