



Reference: PTXBC107129

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EPM7; KV3.1; KV4; NGK2; potassium voltage-gated channel subfamily C member 1; potassium channel, voltage gated Shaw related subfamily C, member 1; voltage-gated potassium channel protein KV3.1; voltage-gated potassium channel subunit Kv4		

For research use only.

GENE DESCRIPTION	potassium voltage-gated channel, Shaw-related subfamily, member 1
GENE ID	3746
ORIGIN SPECIES	Human
ACCESSIONS	BC107129
VECTOR	pDONR223
NCBI SYMBOL	KCNC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr