

GATC-glutamyl-tRNA(Gln) amidotransferase, subunit C homolog (bacterial) Gene

Reference: PTXBC107145

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	15E1.2; glutamyl-tRNA(Gln) amidotransferase subunit C, mitochondrial; glu-AdT subunit C; glutamyl-tRNA(Gln) amidotransferase, subunit C homolog; glutamyl-tRNA amidotransferase subunit C		

For research use only.

GENE DESCRIPTION	glutamyl-tRNA(Gln) amidotransferase, subunit C homolog (bacterial)
GENE ID	283459
ORIGIN SPECIES	Human
ACCESSIONS	BC107145
VECTOR	pDONR223
NCBI SYMBOL	GATC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtcgcggttggtgtgctggccttcgggcccctctgggcccgcgagggccttcacctccaaggcggatcctcaggcagtgccggatcacggctcgggtgatcgagcacctggagcgtctagcgttggacttcggcagccgcgaggcagtgccgcgactggagaa
agctatcgcttcgcccagcgctacgcgcgtggacacagacgggtggagcccatggaatcggtcctggaggacagatgtctatactgagatccgacaatgtggtagaaggcaactgtgctgatgaattactacaaaactcccatcgcgtcgtggaggagtactttgtgcccccc
aggtaatatctcttgcgaagctggatgaacaagagccattccacacagctga