

PAX4-paired box 4 Gene

Reference: PTXBC107149

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KPD; MODY9; paired box protein Pax-4; paired box gene 4; paired domain gene 4; paired box 4		

For research use only.

GENE DESCRIPTION	paired box 4
GENE ID	5078
ORIGIN SPECIES	Human
ACCESSIONS	BC107149
VECTOR	pDONR223
NCBI SYMBOL	PAX4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctgtgctcagtaagctgtgtgttctttagctcctcctgtccaggacactctcttcttgcctatcatgctcactgtcctgtctcactgactcaggtatctaatggctgtgtgagcaagatctagggcgttactaccgcacaggtgtcttggagccaaagggcattgg
gggaagcaagccacggctgtgtacacccctgtgtgtgtcgaattgccagctgaagggtgagtgctcagccctcttgcctgggaatccaacgccagcttgtgtgaagggtttgcacccaggacaagactccaggtgtctcctccatcaaccgagtcctgcggcattacagga
ggaccagggaactaccgtgcacacggctcaggtcaccagctgtttgtgtccagctgtctcactcctccatagtggtctctgagactccccgggtaccacccagggaccgcccaggaatcggaactatcttccccaaagccaagcagagcactggagaagagttccagcgtggg
cagtatctgattcagtgcccggtgaaagctggctactgccactctctgctgaggacacgggtgaggtgtgtgtttccaacagaagagccaatggcgtcggaagagaagctcaagtggaatgcagctgccaggtgtctccaggggtgactgtaccaaggttgcaccagg
aatcatctctgcacaggcccccacagtggttcatgtcaagaaagtactaa