



Reference: PTXBC107157

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PGLYRPIbeta; PGRP-Ibeta; PGRPIB; SBBI67; peptidoglycan recognition protein 4; PGRP-I-beta; peptidoglycan recognition protein I-beta; peptidoglycan recognition protein intermediate beta		

For research use only.

GENE DESCRIPTION	peptidoglycan recognition protein 4
GENE ID	57115
ORIGIN SPECIES	Human
ACCESSIONS	BC107157
VECTOR	pDONR223
NCBI SYMBOL	PGLYRP4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]