

EFNA4-ephrin-A4 Gene

Reference: PTXBC107483

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EFL4; EPLG4; LERK4; ephrin-A4; LERK-4; eph-related receptor tyrosine kinase ligand 4; ligand of eph-related kinase 4; ephrin A4		

For research use only.

GENE DESCRIPTION	ephrin-A4
GENE ID	1945
ORIGIN SPECIES	Human
ACCESSIONS	BC107483
VECTOR	pDONR223
NCBI SYMBOL	EFNA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctgctgccccctgctgaggaactgtctctggcgccgcttctcctggctccccctgctgctggcgccgctcagcctccgacgtactgtggaactccagtaacccccaggttgcttcgagagagcgccgtggtggagctggcgctcaacgattacatgacattgtctgccccact
acgaaggccagggccccctgagggccccgagagcttcttcttgcacatggtggaactggcaggtctatgagtcctccagggcagagggccccggcgctacaagcgtgggtgtgtctccctgccccttggccatgttcaattctcagagaagattcagcgcttcacaccttctccctcgg
cttggagttctacctggagagacttactactacatctcgtgccccactccagagagttctggccaggtctgagggctccaggtgtctgtctgtcgaaggagaggaagctgagtcagccatctgttgggagccctggagagagtgccacatcagggtggcgaggggggacactcc
agccccctctgtcttcttactactgtcttctgattcttctgtcttctgcaattctgtga