

CNPY1-canopy 1 homolog (zebrafish) Gene

Reference: PTXBC107598

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein canopy homolog 1; canopy 1 homolog; canopy FGF signaling regulator 1		
For research use only.			

GENE DESCRIPTION	canopy 1 homolog (zebrafish)
GENE ID	285888
ORIGIN SPECIES	Human
ACCESSIONS	BC107598
VECTOR	pDONR223
NCBI SYMBOL	CNPY1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacgactacaagcttgaggaagaccctgtgacgaaggagagaactttcaagagattcgtctctaggaaaggagacaaaatataccaagaatttaaaaaattgtattttattctgatgcttacagaccttgaaatttgcgtgtgaaactataagaagagtatgaagatgaaatcctcac
ttatcgccaggagacacactatctagctgacaagctgtgcagtgaaaaatcagatctgtgtgaaactctgctaatactagctctag