



Reference: PTXBC107812

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IKK-E; IKK-i; IKKE; IKKI; inhibitor of nuclear factor kappa-B kinase subunit epsilon; I-kappa-B kinase epsilon; IKK-epsilon; IKK-related kinase epsilon; inducible I kappa-B kinase; inducible IkappaB kinase; inhibitor of kappa light polypeptide gene enhanc		

For research use only.

GENE DESCRIPTION	inhibitor of kappa light polypeptide gene enhancer in B-cells, kinase epsilon
GENE ID	9641
ORIGIN SPECIES	Human
ACCESSIONS	BC107812
VECTOR	pDONR223
NCBI SYMBOL	IKBKE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr