

LOC338809-hypothetical protein LOC338809 Gene

Reference: PTXBC107813

| | | | |
|----------------|------|---------------------|----------------|
| Product type: | DNA | Storage Condition: | -20°C |
| Size/Quantity: | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms | | | |

For research use only.

| | |
|-------------------------|---|
| GENE DESCRIPTION | hypothetical protein LOC338809 |
| GENE ID | 338809 |
| ORIGIN SPECIES | Human |
| ACCESSIONS | BC107813 |
| VECTOR | pDONR223 |
| NCBI SYMBOL | LOC338809 |
| LENGTH | ProteoGenix bp |
| SEQUENCE PRIMERS | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

ORF (Only for reference, product provided without STOP codon)

atggagaaaacagagtcattctgtccagaggtgccacccaagactgtggagcctctctcggccctcctgaggagcctgcacaaagaaccagggagctctccagttgaccggcaagccccagggccgcatctccactgccactttgaggagattaaggaccctggcctgtgg
gacaaggcaggatgcttggcaggtgaccacctggggaagctggggagctcctgtggcttccatgttaccttcgaagagcctgccaggaagcccaagattctcacactgtctcacccttgagactccactcaagattgacctgaacctgaaagggcctgaatcagctgac
tactggagcccaaacccggcccactgacaagtatctccctctgagcttcagcctgtcaatgaaggtccctcaccagcctctctcggcagcaagaaggccactctgccagccccaccctacgacacctagctcaggagaagaattgcacccatccaggtgtatgcat
ttattttccggtgctatgacatctgtga