

OR2J2-olfactory receptor, family 2, subfamily J, member 2 Gene

Reference: PTXBC108898

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OR6-19; OR6-8; OR6.3.8; ORL684; dJ80I19.4; hs6M1-6; olfactory receptor 2J2; olfactory receptor 6-8; olfactory receptor OR6-19; olfactory receptor family 2 subfamily J member 2		

For research use only.

GENE DESCRIPTION	olfactory receptor, family 2, subfamily J, member 2
GENE ID	26707
ORIGIN SPECIES	Human
ACCESSIONS	BC108898
VECTOR	pDONR223
NCBI SYMBOL	OR2J2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgattaaaaaatgcaagtcggaagactcttattctacttggattttctaattggcctcagctggaagtagttcttctgtgttattcttctcctgatgacactgacaggaaacctgttcatcatcctgtcctacgtggactcccatctccacacacaaatgtacttctcttca
aacctctcatttctggatctctctacaccaccagctctatccctcagttgtgtggaatctccgggcccggaaaaagaccatctcgtatgctgtgtgcatggttcaactttactttgttctgactgggaatcgagagtgctcctactgggtggtgatgcatatgctgttagctgtgttag
accttgcattacactgtcctcatgcaccctcgttctgccacttgtgtgtggtggttcttgggttaattggtttactatctcagcacttcatctccttacttcttgggtacccttgggacatcgcttagtggaatgacagcacttctcgtttatcatgtgttgacacccatg
caaatgagctgacctcatggtcatgagctccattttgttctacactctcattctcactacatggtgccattgccgggctgtactgagcatgcaatcaaccactgggctcagaaagtgtttaggacatgtggagcccatcttattggttgcctctcttttcatccagtcagtcgcat
gtatccagccacatcagaaattctctgatcagggaagttcattgcccttttatactgtgtcacaccgagtttaactcttaactctacacttcagaaacaagcatgtaaaaggggcagcgaaagagactattgggtgggagtgagggaagtga