

C17orf55-chromosome 17 open reading frame 55 Gene

Reference: PTXBC108933

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf55; long intergenic non-protein coding RNA 482		
For research use only.			

GENE DESCRIPTION	chromosome 17 open reading frame 55
GENE ID	284185
ORIGIN SPECIES	Human
ACCESSIONS	BC108933
VECTOR	pDONR223
NCBI SYMBOL	C17orf55
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgagatccctgaagcgggcccagctccgttcactgctccccggcctcctgtgtttccacacgcagccatgtacagggcagctcatccaccaccccccaaccgcagctccaggggatgtacacctgcctcaggggtgtgcctgggcccctggctgtggaggccacagccgtcaggggaagccctcagctgcaggatgccgtgctcagagaccaaccagccctcgaaggcactgtggccggccagatgtcagcagctccagcaattaggctgggacagatggctcctgtgacacaaggggactgtggggccacaggggaccctgtcacctggaccta ccgtggaggccagggaggagatggactcgaaggcgggaaggccagagaaggcacatttgcagagcagcggccgcatftccagagctcgggagcccaacagggaatcccgtctggccatggggaccaccgcccctggggctgggggatgccgcaggggacggcagagggcagacagggccagggaaaaggaggaggcagagggcagacagggcaggaagagcgccctgcaagtgcccaagggaaggacccaaccagggccctggaccagagcggccgctgtgtggggcaggctgggaaggcccaaggccagcggcgaaggaggagcaggtgagggaccctgtgggcatctgtgggaacaaggccacgtttctccctgcgcagcgtttaatacaggacattcctgtgtgttccaaaagtgtgcacactacatgggtgtcgtga