

TM2D2-TM2 domain containing 2 Gene

Reference: PTXBC109050

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BLP1; TM2 domain-containing protein 2; BBP-like protein 1; beta-amyloid-binding protein-like protein 1; TM2 domain containing 2		

For research use only.

GENE DESCRIPTION	TM2 domain containing 2
GENE ID	83877
ORIGIN SPECIES	Human
ACCESSIONS	BC109050
VECTOR	pDONR223
NCBI SYMBOL	TM2D2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgctaggtggtgcccgttagtacttctgtcgccgagcgcgcttctgctggtggaatttacttctgctgcaattgtgtgtctcgagccactgcacaaatgcgaccgctgagcctgagctcacatccgctggcgccgcccagccggaggcccccgggggtgctgcgagctg
ggaatatggcgaccgccactctccgggtcatcctctgtcttacctacctgafgaatttatagaatgtgaagaccagtgatcatgttggaatgcaactgcatccagggaacttggttatggtgtctcaagtgcggtcagcctacagcgacgtggaaacacacttcagtcagtcagtcag
ccttagatggaattgagtggtcctaggaaccttctacgagaaaataaacctgtataaagtataccggacactacttcataaacactttactctactctcttctcgggatgttttggtgtggaatcattctgtttgggacacactggcactgcagtagggaagctgttgacgttgaggac
ttgggatttggtgtgttgaccttatttgcataattactggaggctgatccaagtgatggcagcaactggtgcactgtttactaa