



Reference: PTXBC109126

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCDC147; C10orf80; bA127L20.4; bA127L20.5; bA554P13.1; cilia- and flagella-associated protein 58; bA127L20.4 (novel protein); bA127L20.5 (novel protein); coiled-coil domain containing 147; coiled-coil domain-containing protein 147; cilia and flagella asso		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 147
GENE ID	159686
ORIGIN SPECIES	Human
ACCESSIONS	BC109126
VECTOR	pDONR223
NCBI SYMBOL	CCDC147
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]