



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MAGEA5-melanoma antigen family A, 5 Gene

Reference: PTXBC109187

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT1.5; MAGE5; MAGEA4; melanoma-associated antigen 5; MAGE-5 antigen; MAGE-5a antigen; MAGE-5b antigen; cancer/testis antigen 1.5; cancer/testis antigen family 1, member 5; melanoma antigen family A, 5; melanoma antigen family A5; MAGE family member A5		

For research use only.

GENE DESCRIPTION	melanoma antigen family A, 5
GENE ID	4104
ORIGIN SPECIES	Human
ACCESSIONS	BC109187
VECTOR	pDONR223
NCBI SYMBOL	MAGEA5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctctgagcagaagagtcgacgactgcaagcctgagggaagccttgacaccaagaagaggccctggcgctgggtgtgcaggctgccactactgaggagcaggagctgtgtctctctctctctctggtccaggacccttgggggaggtgctgctgtggttcaccg
ctctctcaagagtcttcaggaggagctctcgcatccccactgccaatgatttactctatggaggcaatcattaaagggctccagcaaccaagaagaggaggggccaagcactccctgaccagagctgtgtgtccgagcagcactcagtaagaaggtggtgacttgatttcttgc
gttcctaaagtattaa