



Reference: PTXBC109209

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	P2RXL1; P2X6; P2XM; P2X purinoceptor 6; ATP receptor; purinergic receptor P2X, ligand gated ion channel, 6; purinergic receptor P2X-like 1, orphan receptor; purinoceptor P2X6; purinoreceptor P2X6; skeletal muscle-expressed P2X; purinergic receptor P2X 6		

GENE DESCRIPTION	purinergic receptor P2X, ligand-gated ion channel, 6
GENE ID	9127
ORIGIN SPECIES	Human
ACCESSIONS	BC109209
VECTOR	pDONR223
NCBI SYMBOL	P2RX6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1