



Reference: PTXBC109210

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	THETA; gamma-aminobutyric acid receptor subunit theta; GABA(A) receptor subunit theta; gamma-aminobutyric acid (GABA) A receptor, theta; gamma-aminobutyric acid type A receptor theta subunit		

For research use only.

GENE DESCRIPTION	gamma-aminobutyric acid (GABA) receptor, theta
GENE ID	55879
ORIGIN SPECIES	Human
ACCESSIONS	BC109210
VECTOR	pDONR223
NCBI SYMBOL	GABRQ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]