

STAC2-SH3 and cysteine rich domain 2 Gene

Reference: PTXBC109231

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	24b2/STAC2; 24b2; SH3 and cysteine-rich domain-containing protein 2; SRC homology 3 and cysteine-rich domain-containing protein 2; SH3 and cysteine rich domain 2		

For research use only.

GENE DESCRIPTION	SH3 and cysteine rich domain 2
GENE ID	342667
ORIGIN SPECIES	Human
ACCESSIONS	BC109231
VECTOR	pDONR223
NCBI SYMBOL	STAC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcaaaagtcagctccacctctgtgtctctgaggagatctccaccagcaatgccaggcaagacgtccacctcttccgcccaacttcagttccctctctctgtgtcatgagccgccaccagctctgtccacaagcaaaagagtcacccaccactggggacagtgggaaggtggac
ccgtgtctacgagacctgcgtatggcactccctggcactgatgaaccgtccagtttcagcagcacctctgagtcctccgacaaggagcctgagtgagcgggaatgagctgaccgaggatggggaaggcagcaltccgagctctgaggaggccctgtgtgacagtgcattccagttatt
cacagcccccagcagagagtgaaggccaggaccagaggagaagagtcctggacagcagctcccaagccacctgcggaaggatgtggggccatgtactcctacgttgacacttacaagtcttctgccccaggagaacaatgatctgtgctctgcagcctggagatcggtatcatgctg
gtggatgacttaacgaggactgtgtggaaggccaagatcggcgaccgggtgtgtcttccagctaatgtgtgcaacgggtgagggccaggcgagaatgtttggcgctgtgccaaccttctccgggaacaaggacaagggttacatgagcctcaaggagaaccagatctgcgtggg
cggtgggcagaagcaaggatgctgacggcttcatccgctgagcagtggtgcaagaagcggggcctgtgtgccagtcgacgacctgactgagatctga