

OR2B6-olfactory receptor, family 2, subfamily B, member 6 Gene

Reference: PTXBC109251

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OR2B1; OR2B1P; OR2B5; OR2B6P; OR5-40; OR5-41; OR6-31; dJ408B20.2; olfactory receptor 2B6; hs6M1-32; olfactory receptor 2B1; olfactory receptor 2B5; olfactory receptor 5-40; olfactory receptor 6-31; olfactory receptor OR6-4; olfactory receptor, family 2, s		

For research use only.

GENE DESCRIPTION	olfactory receptor, family 2, subfamily B, member 6
GENE ID	26212
ORIGIN SPECIES	Human
ACCESSIONS	BC109251
VECTOR	pDONR223
NCBI SYMBOL	OR2B6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaattgggtaaatgacagcatcacaggagttattctgctgggtttctcagatcgacctggctggagttccactcctgtggtctcttgattcttacactgtgacctctttggcaatctgaccattattctagtgtcacgctggacaccaaactcatacccccatgtattttttctaccaat
ctatcactcctggatctttgttacaccatgtacagtcaccacaatgctagtaaatgtatgcagcatcaggaaagtaacagttatcgtggctgtgtgagccagctttcatatttctggccttgggggctactgaatatcttctctggcgtcatgtcctttgataggtttgtagctattgtcgccctc
tccattactcagttatcatgcaccagagactctgctccagttggcagctgcatcctgggttactggttttagtaactcagtggtgtgtctaccctgactctccagctgcccactctgtgacccctatgtgatagatcactttctctgtgaagtcctgcaactgctcaagttatctgtgttgagacaaca
gcaaatgaggctgaactattcctgtcagtgagcttccatcaataaccctgacactatccttatatcatatgcttttattgtccgagcagttaggagatacagctctgctgaaggctgacaaaaagcatttgggacatgtgtgtcccatctaatgtgtgtctctttttatagtagcagcgtctctg
tgtacctgcaaccaccttgcgccagctccaaggaccaaggaaagatggtttctctctctatggaatcattgacccatgctgaatcccttatatatatacttaggaacaaggaggttaaaagggttgggtgcaagagctcttcaatcaagaataaa