



Reference: PTXBC109297

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CALDAG-GEFI; CALDAG-GEFII; RASGRP; hRasGRP1; RAS guanyl-releasing protein 1; RAS guanyl nucleotide-releasing protein 1; RAS guanyl releasing protein 1 (calcium and DAG-regulated); calcium and DAG-regulated guanine nucleotide exchange factor II; calcium- a		

For research use only.

GENE DESCRIPTION	RAS guanyl releasing protein 1 (calcium and DAG-regulated)
GENE ID	10125
ORIGIN SPECIES	Human
ACCESSIONS	BC109297
VECTOR	pDONR223
NCBI SYMBOL	RASGRP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]