



Reference: PTXBC109308

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMEM16D; anoctamin-4; transmembrane protein 16D; anoctamin 4		
For research use only.			

GENE DESCRIPTION	anoctamin 4
GENE ID	121601
ORIGIN SPECIES	Human
ACCESSIONS	BC109308
VECTOR	pDONR223
NCBI SYMBOL	ANO4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr