



Reference: PTXBC109311

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	rab effector MyRIP; SLAC2-C; SLAC2C; Slp homologue lacking C2 domains; exophilin-8; myosin-VIIa- and Rab-interacting protein; slp homolog lacking C2 domains c; synaptotagmin-like protein homologue lacking C2 domains-c; synaptotagmin-like protein lacking C		

GENE DESCRIPTION	myosin VIIA and Rab interacting protein
GENE ID	25924
ORIGIN SPECIES	Human
ACCESSIONS	BC109311
VECTOR	pDONR223
NCBI SYMBOL	MYRIP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]