

SPINK7-serine peptidase inhibitor, Kazal type 7 (putative) Gene

Reference: PTXBC109385

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ECG2; serine protease inhibitor Kazal-type 7; ECRG-2; esophagus cancer-related gene 2 protein; esophagus cancer-related gene-2; serine peptidase inhibitor, Kazal type 7 (putative)		

For research use only.

GENE DESCRIPTION	serine peptidase inhibitor, Kazal type 7 (putative)
GENE ID	84651
ORIGIN SPECIES	Human
ACCESSIONS	BC109385
VECTOR	pDONR223
NCBI SYMBOL	SPINK7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagatcactgggggtctccttctgtctgtacagtgggtctattctgtagcagctcagaagctgctagctgtctccaaaaaagtggactgcagcattacaagaagtatccagtgggtggccatcccctgccccatcacatactaccagttgtgggtctgactacatcacctatgggaatg
aafgtcacttgtgaccgagagcttgaagaaagtaatggaagagttcagtttcttcacgatggaagttgctaa