

C14orf142-chromosome 14 open reading frame 142 Gene

Reference: PTXBC109390

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C14orf142; C14orf142; PNAS-127; GON7, KEOPS complex subunit homolog		

For research use only.

GENE DESCRIPTION	chromosome 14 open reading frame 142
GENE ID	84520
ORIGIN SPECIES	Human
ACCESSIONS	BC109390
VECTOR	pDONR223
NCBI SYMBOL	C14orf142
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacattaaaagatggctgttgccttattgtcacactaggacatttaagggaagacttatgctctaatttggaaaaagcatttgtgtctcctgggacaatttgcacaagagaaaactggttttctgttctagcaaagatagtaaaggcttagaataaaatctgtttcttgagaaataaattcagaccctggaagaaaatag