

**SPINK4-serine peptidase inhibitor, Kazal type 4 Gene**

Reference: PTXBC110068

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL136; PEC-60; PEC60; serine protease inhibitor Kazal-type 4; epididymis luminal protein 136; gastrointestinal peptide; peptide PEC-60 homolog; serine peptidase inhibitor, Kazal type 4		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	serine peptidase inhibitor, Kazal type 4
<b>GENE ID</b>	27290
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC110068
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	SPINK4
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atggcgtccgccagtgggtaatgccttggccttggctgccctctgttggagcagggaagtgccagtggcagcaggaaagtccttctcaagaatgccatcttgaacacatggtagatgtccaacctgttccagatgtccaacctgggtctgggacctgatggctcacat  
acgaatgaatgccagctcttggccggataaaaaccaaacaggacatccagatcatgaaagatggcaaatgctga