

RAB11B-RAB11B, member RAS oncogene family Gene

Reference: PTXBC110081

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB11B, member RAS oncogene family; RAB11B, member of RAS oncogene family; H-YPT3; ras-related protein Rab-11B; GTP-binding protein YPT3		

For research use only.

GENE DESCRIPTION	RAB11B, member RAS oncogene family
GENE ID	9230
ORIGIN SPECIES	Human
ACCESSIONS	BC110081
VECTOR	pDONR223
NCBI SYMBOL	RAB11B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggacccgggacgacgagtagctacattcaaaagtgtgtctatcggggactcaggcgtgggcaagagcaacctgtgtcgtcgttccccgcaacgagttcaacctggagagcaagagcaccatcggcgtggagttcggcaccgcagcatccaggtggacggcaagacatcaaggcgcagatctgggacaccgtggccaggagcgtacccgcctacacctccgctactaccgtgtgtcagttggcgccctgtgtgtgtacgacatcgccaagcactgacctatgagaacgtggagcgtgtgctgaaggagctgcgggaccacgcagacagcaacatcgtcatcatgtgtgtgggcaacaagagtgacctgcgccacctgcgggctgtgtccactgacgaggcccgcccttcgagaaaaaacaactgtctctcatcgagacctcagccttggattccactaactgagagggaagcattcaagaacatcctcacagatctaccgcatcgtgtcaca gaaacagatcgagaccgtgtgtccacgacgagtgccccggggaacaacgtgtgtgacatcagcgtgccgccaccacgacggacagaaagcccaacagctgcagtgctgccagaacctgtga