



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

EFNA3-ephrin-A3 Gene

Reference: PTXBC110406

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EFL2; EPLG3; Ehk1-L; LERK3; ephrin-A3; EFL-2; EHK1 ligand; LERK-3; eph-related receptor tyrosine kinase ligand 3; ligand of eph-related kinase 3; ephrin A3		

For research use only.

GENE DESCRIPTION	ephrin-A3
GENE ID	1944
ORIGIN SPECIES	Human
ACCESSIONS	BC110406
VECTOR	pDONR223
NCBI SYMBOL	EFNA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgctgcctgctgctgctgctgctgctgcccgtgcgcgtgctgcgcgtgctgctggcccaagggccggaggggcgctgggaaacacgcaltgcgggtgactggaacagctccaaccagcacctcggcgagagggtctacaccgtgcagctgaacgtgaacgactatc
ggataatttactgcccgactacaacagctcgggggtggggcccgggggcgggaccggggccggaggcggggcagagcagtagctgctgtatcgtgtagccgaacggctaccgcactgcaacgccagcgagggtctcaagcgtgggagtgaacggcgccgacgccccgc
acagcccacataagattctcggagaagtctccagcgctacagcgcccttctctctgggtacagagttccagcgccgcacgagtagtactactacatccacgcccaactacaacctgcactggaagtgtctgaggatgaaggtgttcgtctgctgcgtccacatcgactccgggagaagcc
gtgccacactccccagttaccatggggcccaatgtagaatacaacgtgctggaagactttgaggagagaacccctcaggtgcccaagctgagaagagcatcagcgggaccagcccccacgggaacacgtgccctggccgtgggcaltgccttctctatgactgttgc
ctctagc