

COX19-COX19 cytochrome c oxidase assembly homolog (S. cerevisiae) Gene

Reference: PTXBC110420

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COX19, cytochrome c oxidase assembly factor; COX19 cytochrome c oxidase assembly homolog; cytochrome c oxidase assembly protein COX19; cytochrome c oxidase assembly homolog 19; hCOX19		

For research use only.

GENE DESCRIPTION	COX19 cytochrome c oxidase assembly homolog (S. cerevisiae)
GENE ID	90639
ORIGIN SPECIES	Human
ACCESSIONS	BC110420
VECTOR	pDONR223
NCBI SYMBOL	COX19
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgaccgccatgaatttcggaccagagctccagccgcccgggacagggcagctcccgctggatcacttagtgatgtaaaagctttaaagagaaatcatgaagtgtcttcatacaataatttgaaaatgctttgtgcagaaaggaaatcaaaagaatatttagaatgca
ggatggagagaaaattgatgctacaagaaccattggagaactgggatttggagactgactagtggaaatcagaggcaaaaaatga