



Reference: PTXBC110430

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRSP7; CRSP70; mediator of RNA polymerase II transcription subunit 26; ARC70; CRSP complex subunit 7; activator-recruited cofactor 70 kDa component; cofactor required for Sp1 transcriptional activation subunit 7; cofactor required for Sp1 transcriptional		

For research use only.

GENE DESCRIPTION	mediator complex subunit 26
GENE ID	9441
ORIGIN SPECIES	Human
ACCESSIONS	BC110430
VECTOR	pDONR223
NCBI SYMBOL	MED26
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]