



Reference: PTXBC110445

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HUT2; UT-A2; UT2; UTA; UTR; hUT-A6; urea transporter 2; solute carrier family 14 (urea transporter), member 2; urea transporter, kidney; solute carrier family 14 member 2		

For research use only.

GENE DESCRIPTION	solute carrier family 14 (urea transporter), member 2
GENE ID	8170
ORIGIN SPECIES	Human
ACCESSIONS	BC110445
VECTOR	pDONR223
NCBI SYMBOL	SLC14A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]