

ZNF137-zinc finger protein 137 Gene

Reference: PTXBC110455

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	zinc finger protein 137
GENE ID	7696
ORIGIN SPECIES	Human
ACCESSIONS	BC110455
VECTOR	pDONR223
NCBI SYMBOL	ZNF137
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgtggcaaggttttagtagaaaagcacacttgcattgcatcatagactcactgtctaaaggttttcaatcaaatcaaaccttgcacacatcagagagtttactggagagaaaccttacaagtgaatgagtggtggcaaaagccttaagtgggaagtcgtcactttttatcatcaa
gcaatccatggtgtagggaaactttgcaaatgtaattgtcacaaagcttcagtaattgctacaaccttgcacacactggagaatccgtaattgaagacagatcttacaagtgttaataaatgtgtaaaattttcagacatcaatcatacttgcagtttatcagcgaactcactggagagaa
accttacaatatcatgactgtggcaaggtcttcagttcaagcttcctcatgcaaaacataggagaattcatcagagagaaacctcacaagtgtgatgattgtggcaaaagcttgcattcagttcacacctcattagacatcagagaatccatactggacataaatcttacaagtcttaagt
gtggcaaggtcttcagttgtggcactccatgcagaacatcagaaaattcattttga