



Reference: PTXBC110459

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf49; OBTP; TOMM6; prickle-like protein 4; over-expressed breast tumor protein; overexpressed breast tumor protein; overexpressed breast tumour protein; prickle homolog 4; translocase of outer mitochondrial membrane 6 homolog; prickle planar cell polar		

GENE DESCRIPTION	prickle homolog 4 (Drosophila)
GENE ID	29964
ORIGIN SPECIES	Human
ACCESSIONS	BC110459
VECTOR	pDONR223
NCBI SYMBOL	PRICKLE4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]