



Reference: PTXBC110506

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LMBR1 domain-containing protein 2; LMBR1 domain containing 2		
For research use only.			

GENE DESCRIPTION	LMBR1 domain containing 2
GENE ID	92255
ORIGIN SPECIES	Human
ACCESSIONS	BC110506
VECTOR	pDONR223
NCBI SYMBOL	LMBRD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr