



Reference: PTXBC110518

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein GPR107; GCDRP; LUSTR1; bA138E2.2; lung seven transmembrane receptor 1; G protein-coupled receptor 107		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 107
GENE ID	57720
ORIGIN SPECIES	Human
ACCESSIONS	BC110518
VECTOR	pDONR223
NCBI SYMBOL	GPR107
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr